

Montgomery County Fire and Rescue Service Division of Operations Emergency Medical and Integrated Healthcare Services

Office of Medical Oversight Clinical Practice Guideline

A. Initial Ventilator Settings Quick-Reference

Patients 13 years and older

Mode: Volume AC

Tidal Volume: 7 mL/kg

• Respiratory Rate: 16 breaths/min

• PEEP: 8 cm H2O

FiO2: 0.5

Inspiratory Time: 1 secondInspiratory Trigger: -5 cm H2O

Patients aged 9-12 years

Mode: Volume AC

• Tidal Volume: 6 mL/kg,

• Respiratory Rate: 20 breaths/min

PEEP: 5 cm H2O

• FiO2: 0.5

Inspiratory Time: 0.75 seconds
Inspiratory Trigger: -5 cm H2O