

B. Blood Administration Bag Inventory

EMSDO Vehicle

1. DELTA ICE Cooler

- PCM Pack
- Low Titer O+ Whole Blood unit(s)
- Portable Power Bank with Cable

2. Blood Administration Bag

- QinFlow Warrior Lite Blood Warmer
- QinFlow Warrior Lite Battery
- QinFlow Warrior Lite Blood Warmer Extension Cable
- QinFlow Warrior Lite Blood Warmer Compact Disposable Units (CDU) (2)
- 250ml Normal Saline (2)
- LifeFlow Infuser with Filtered “Y” Blood Tubing
- EZ-IO drill with IO needles, adapter, alcohol preps, and securing device
- 2.5-gallon bag
- Black sharpie
- Laminated instruction pages
- Patient Bracelets (2)
- Trauma shears
- Carabiners (2)

Deployment Site

- Freezer, Set to -16°C
- PCM Packs inside the freezer (2)
- Moisture Collection Mat
- Towel to dry PCM Pack
- QinFlow Warrior Lite Blood Warmer Batteries (2)
- Portable Power Banks (2)